MASSAGE GUN

U head
for vertebrae

Small round head
for joints

Ball head
for large muscle groups

Small flat head
for all parts of the body

EN INSTRUCTIONS FOR USE - MASSAGE GUN

Instructions

1. On/off button (press and hold to start / press and hold to turn off)
2. Press the power button to change gears
3. Removable and replaceable vibrator. When disassembling the vibrator, just pull it out with force.
4. Please choose a 5V charger for charging
5. For the first time, please charge it for 3 hours. To charge, please connect the data cable to the 5V adapter.
6. Button type: the 6 lights flash during charging. When fully charged, all 6 lights light up at the same time. During operation, the lights indicate the gear.
7. Display type: when charging, the display panel shows the electric charge, during use, the gear.

Safety instructions

Hazard:

To reduce the risk of electric shock, fire and personal injury, this product must be used in accordance with the following instructions.

For adults only

It can only be used on dry, clean skin surfaces of the body, it cannot be used through clothes. Gently press and move on the skin. The duration for each body part is as described below.

Use the device only on the soft tissues of the body without causing pain or discomfort.

Do not use on the head or any hard or bony part of the body.

Use only the vibrating head that can produce the best results.

It can cause bruising, either under moderate or normal pressure.

Please always check the area receiving physiotherapy.

Discontinue use immediately if any signs of pain or discomfort occur.

Keep fingers, hair or other body parts away from the shaft and back of the vibrating head as they may become pinched.

Do not insert any object into the ventilation opening of the device.

Do not submerge it in water or allow water to enter the ventilation opening of the device.

Do not drop or misuse the device.

Only charge the device with the appropriate direct voltage (DC)..

Carefully check the integrity of the device and battery before each use.

Do not alter or modify the device in any way.

Never operate or charge the device without supervision.

Indications and contraindications

By massaging the deeper muscle fibers, you can reduce the pain resulting from joint sprains, bursitis, muscle spasms, neuralgia and other musculoskeletal disorders in a short time. Almost all soft tissue disorders are related to trauma, muscle overuse, or abnormal posture, and most of the scar tissue involved in these disorders can be treated with this tool. When using the device, avoid excessive strain and prolonged physiotherapy. Also avoid sensitive areas such as the head, face and neck vertebrae. Discontinue physical therapy if inflammation, swelling, or increased pain occurs. Ongoing research is about the effect of vibration training on certain health conditions. This is likely to reduce the types of contraindications below, and practical experience shows that in many cases the physiotherapy recommendation includes the inclusion of vibration training in physiotherapy programs. This must be done by or in the presence of a doctor, specialist or physiotherapist.

Indications

- Pain and spasms due to muscle injury, sprain, strain

- It helps the flow of edema fluid in the swollen area
- Loosening of thickened connective tissues and fascia
- It reduces the accumulation of lactic acid in the muscles
- Increases joint mobility
- Eliminates muscle fatigue

Contraindications

- Aneurysm, bleeding, use of blood thinners.
- Heart disease, use of a pacemaker or defibrillator
- Pregnancy, cancer
- The device cannot be used in the event of a bone fracture or fixation surgery within 90 days
- Sensitive area: head, face, cervical vertebrae, spine
- Areas with less muscle are positions close to bones: tibia, back of foot, back of hand

Method of use

Functional program of the lower limbs

Supine position

Small flat head - vertically down

Knee ligament

Biceps femoris (Thigh flexor muscle)

20 sec-30 sec x 3 times

Semitendinosus (Posterior hamstring)

20 sec-30 sec x 3 times

Semimembranosus (horse muscle)

20 sec-30 sec x 3 times

Direction of movement:

from top to bottom.

Small flat head - outward - 45°

Gracilis (inner leg)

10 sec-15 sec x 2 times

Adductor longus

10 sec-15 sec x 2 times

Adductor brevis (short thigh adductor muscle)

10 sec-15 sec x 2 times

Triceps surae

Gastrocnemius (feet)

20 sec-30 sec x 3 times

Soleus (leg flexor muscle)

20 sec-30 sec x 3 times

Direction of movement:

from top to bottom

Small round head - vertically down

Popliteus (knee bend)

Knee straightening position

15 sec-30 sec

Knee bending position

15-30 sec

Method of use

Knee joint

Supine position

Small flat head - vertically down

Quadriceps femoris

Rectus femoris (straight thigh muscle)

30 sec-45 sec x 2 times

Vastus lateralis (external muscle)

30 sec-45 sec x 2 times

Vastus medialis (inner vasculature)

30 sec-45 sec x 2 times

Vastus intermedius

30 sec-45 sec x 2 times

Direction of movement:

from top to bottom

The patient slowly and slightly bends the knee.

Small round head - vertically down

The edge of the hip bone - the front part of the hip bone

30 sec-1 min

Rectus femoris, inguen (rectus femoris, hip)

1 minute

Tensor fasciae latae

15 sec-20 sec x 2 times

Patellar ligament area (patella tendon)

20 seconds x 3 times

Supine position

Small round head - vertically down

Knee bend

30 sec-45 sec

Knee straightening position

15 sec-30 sec

Knee bending position

15 sec-30 sec

Triceps surae

Gastrocnemius (feet)

20 sec-30 sec x 3 times

Soleus (leg flexor muscle)

20 sec-30 sec x 3 times

Method of use

Plantar fasciitis Supine lying (one sided)

Small flat head - vertically down

Knee ligament

Biceps femoris (thigh flexor muscle)

20 sec-30 sec x 3 times

Semitendinosus (hamstring)

20 sec-30 sec x 3 times

Semimembranosus (horse muscle)

20 sec-30 sec x 3 times

Direction of movement:

from top to bottom.

U head vertically down

Knee bend

Knee straightening position

15 sec-30 sec

Knee bending position

15 sec-30 sec

Triceps surae

Gastrocnemius (feet)

20 sec-30 sec x 3 times

Soleus (leg flexor muscle)

20 sec-30 sec x 3 times

The sole of the front foot touches the ground

Small round head - vertical to the planar pedis

Plantar fascia

Functional program of the sacroiliac part

Supine (one-sided)

Small flat head - vertically down

Gluteus medius

30 sec-1 minute x 4 times

Gluteus maximus (large gluteal muscle)

30 sec-1 minute x 4 times

Piriformis (pear-shaped muscle)

1-1.5 minutes

Articulatio sacroiliaca (sacral joint)

1-1.5 minutes

Direction of movement:

from the inside up, from the outside down, along the direction of the muscle

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